

# Quick Start Guide



1. Download and install Yunmai app from app or play store or <http://www.iyunmai.com/en/app.html> , login and setup your profile with your information (name, gender, height, age).
2. Place the scale on a hard flat surface like tiles and clean the scale surface with a wet cloth. Check that the scale surface and your feet are not too dry as this can have an effect on the body fat reading.
3. Turn on the Bluetooth on your cellphone, open and login Yunmai app and step on the scale with barefoot, the scale will connect to your phone automatically without a password needed, you might need manually pair it with passkey "000000". You may need to turn on location service and restart your phone to connect to the scale. Please reinstall or change new batteries properly if cannot connect to Bluetooth successfully. Please also update your iOS or Android to latest version and reinstall the latest Yunmai app.
4. Do not move until the scale displays weight and Bluetooth logo "B". Yunmai app will show "connected" on the top instead of "connecting..." and the Bluetooth logo light will turn on, indicating the Bluetooth connection is successful while the app shows all body measurements.
5. Scale might need calibrated to show accurate readings. Step on and step off the scale quickly, the scale will show 0.0 lbs, which means calibration is complete.
6. Please make sure the previous device is disconnected (just close the Yunmai app) before pairing another phone with the scale. The scale can only be connected to one device at a time. If you have several Yunmai scales, please remove batteries from the one you do not want to use to disconnect it with the phone.
7. If you have multiple users sharing one app, make sure tap the profile picture on the top left and switch to the correct member before stepping on the scale.